



Product Name: Recombinant Human GAS6 Protein
Catalog #: PEH5486



Summary

Name	GAS6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q14393
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FLRPRQRRAF QVFEEAKQGH LERECVEELC SREEAREVFE NDPETDYFYP RYLDCINKYG SPYTKNSGFA TCVQNLDPQC TPNPCDRKGT QACQDLMGNF FCLCKAGWGG RLC DKDVNEC SQENGGCLQI CHNKPGSFHC SCHSGFELSS DGRTCQDIDE CADSEACGEA RCKNLPGSYS CLCDEGFAYS SQEKACRDVD ECLQGRCEQV CVNSPGSYTC HCDGRGGLKL

image



Background

Alternative Names Growth arrest-specific protein 6 (GAS-6) (AXL receptor tyrosine kinase ligand)

Note

Product Name: Recombinant Human GAS6 Protein
Catalog #: PEH5486



For Research Use Only , Not for Diagnostic Use