


Product Name: Recombinant Human PRSS2 Protein
Catalog #: PEH5490



Summary

Name	TRY2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P07478
Host	E.coli
Species	Human
Molecular Mass	36 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVGGYICEEN SVPYQVSLNS GYHFCGGSLI SEQWVVSAGH CYKSRIQVRL GEHNIEVLEG NEQFINAAKI IRHPKYNSRT LDNDILLIKL SSPAVINSRV SAISLPTAPP AAGTESLISG WGNTLSSGAD YPDELQCLDA PVLSQAECEA SYPGKITNNM FCVGFLEGGK DSCQGDSGGP VVSNLQGI VSWGYGCAQK NRPGVYTKVY NYVDWIKDTI AANS

image



Background

Alternative Names Trypsin-2 (EC 3.4.21.4) (Anionic trypsinogen) (Serine protease 2) (Trypsin II)

Note

Product Name: Recombinant Human PRSS2 Protein
Catalog #: PEH5490



For Research Use Only , Not for Diagnostic Use