


Product Name: Recombinant Human CADM3 Protein
Catalog #: PEH5492



Summary

Name	CADM3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8N126
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QKPIITGYKS SLREKDTATL NCQSSGSKPA ARLTWRKGDQ ELHGEPTRIQ EDPNGKTFTV SSSVTFQVTR EDDGASIVCS VNHESLKGAD

image



Background

Alternative Names	Cell adhesion molecule 3 (Brain immunoglobulin receptor) (Immunoglobulin superfamily member 4B) (IgSF4B) (Nectin-like protein 1) (NECL-1) (Synaptic cell adhesion molecule 3) (SynCAM3) (TSLC1-like protein 1) (TSLL1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use