

Product Name: Recombinant Human UNC5C Protein
Catalog #: PEH5493



Summary

Name	UNC5C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O95185
Host	E.coli
Species	Human
Molecular Mass	39 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ETFPSDPPEP LPHFLIEPEE AYIVKNKPVN LYCKASPATQ IYFKCNSEWV HQQDHIVDER VDETSGLIVR EVSIEISRQQ VEELFGPEDY WCQCVAWSSA GTTKSRKAYV RIAYLRKTFE QEPLGKEVSL EQEVLLQCRP PEGIPVAEVE WLKNEDIIDP VEDRNFYITI DHNLIQKAR LSDTANYTCV AKNIVAKRKS TTATVIVYVN GGWSTWTEWS VCNSRCGRGY QKRTRCTNP APLNGGAFCE GQSVQKIACT TPCVDGRWT PWSKWSTCGT ECTHWRRREC TAPAPKNGGK DCDGLVLQSK NCTDGLCMQT APDSDDVALY

image



Background

Alternative Names Netrin receptor UNC5C (Protein unc-5 homolog 3) (Protein unc-5 homolog C)

Product Name: Recombinant Human UNC5C Protein
Catalog #: PEH5493



Note

For Research Use Only , Not for Diagnostic Use