


Product Name: Recombinant Human GFRA3 Protein
Catalog #: PEH5496



Summary

Name	GFRA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O60609
Host	E.coli
Species	Human
Molecular Mass	35 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MNSCLQARRK CQADPTCSAA YHHLDSCTSS ISTPLPSEEP SVPADCLEAA QQLRNSSLIG CMCHRRMKNQ VACLDIYWTV HRARSLGNYE LDVSPYEDTV TSKPWKMNLS KLNMLKPDS LCLKFAMLCT LNDKCDRLRK AYGEACSGPH CQRHVCLRQL LTFEKAAP HAQGLLLCPC APNDRGCGER

image



Background

Alternative Names GDNF family receptor alpha-3 (GDNF receptor alpha-3) (GDNFR-alpha-3) (GFR-alpha-3)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human GFRA3 Protein
Catalog #: PEH5496

