


Product Name: Recombinant Human CSF2RA Protein
Catalog #: PEH5499



Summary

Name	CSF2R_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P15509
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	REIRCPYYIQ DSGTHVGCHL DNLSGLTSRN YFLVNGTSRE IGIQFFDSSL DTKKIERFNP PSNVTVRCNT THCLVRWKQP RTYQKLSYLD FQYQLDVHRK NTQPGTENLL INVSGDLENR YNFPSSSEPPRA KHSVKIRAAD VRILNWSSWS EAIEFGSDDG

image



Background

Alternative Names	Granulocyte-macrophage colony-stimulating factor receptor subunit alpha (GM-CSF-R-alpha) (GMCSFR-alpha) (GMR-alpha) (CDw116) (CD antigen CD116)
--------------------------	---

Note

Product Name: Recombinant Human CSF2RA Protein
Catalog #: PEH5499



For Research Use Only , Not for Diagnostic Use