


Product Name: Recombinant Human MDGA1 Protein
Catalog #: PEH5500



Summary

Name	MDGA1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8NFP4
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IKESENIQLG QDLKLSCHVD AVPQEKVTYQ WFKNGKPARM SKRLLVTRND PELPAVTSSL ELIDLHFSYD GTYLCMASFP GAPVPDLSVE VNISSETVPP TISVPKGRAV VTVREGSPAELQCEVRGKPR PPVLWSRVDK EAALLPSGLP LEETPDGKLR LERSRDMMSG TYRCQTARYN GFNVPRPRAQ

image



Background

Alternative Names	MAM domain-containing glycosylphosphatidylinositol anchor protein 1 (GPI and MAM protein) (GPIM) (Glycosylphosphatidylinositol-MAM) (MAM domain-containing protein 3)
--------------------------	---

Note

Product Name: Recombinant Human MDGA1 Protein
Catalog #: PEH5500



For Research Use Only , Not for Diagnostic Use