


**Product Name: Recombinant Human WFIKKN1 Protein**  
**Catalog #: PEH5504**



## Summary

<b>Name</b>	WFKN1_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	N-SUMO & C-His
<b>Accession #</b>	Q96NZ8
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	33 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	LPAVQGPCRG WEPRWAYSPL LQQCHPFVYG GCEGNGNNFH SRESCEDACP VPRTPPCRAC RLRSKLALSL CRSDFAIVGR LTEVLEEPEA AGGIARVALE DVLKDDKMGL KFLGTKYLEV TLSGMDWACP CPNMTAGDGP LVIMGEVRDG VAVLDAGSYV RAASEKRVKK ILELLEKQAC

## image



## Background

<b>Alternative Names</b>	WAP, Kazal, immunoglobulin, Kunitz and NTR domain-containing protein 1 (Growth and differentiation factor-associated serum protein 2) (GASP-2) (hGASP-2) (WAP, follistatin, immunoglobulin, Kunitz and NTR domain-containing protein)
--------------------------	---

**Product Name: Recombinant Human WFIKKN1 Protein**  
**Catalog #: PEH5504**



---

**Note**

For Research Use Only , Not for Diagnostic Use