


Product Name: Recombinant Human IL22RA1 Protein
Catalog #: PEH5511



Summary

Name	I22R1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q8N6P7
Host	E.coli
Species	Human
Molecular Mass	37 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PSDLLQHVKF QSSNFENILT WDSGPEGTPD TVYSIEYKTY GERDWVAKKG CQRITRKSCN LTVETGNLTE LYYARVTAVS AGGRSATKMT DRFSSLQHTT LKPPDVTCIS KVRSIQMIVH PTPTPIRAGD GHRLTLEDIF HDLFYHLELQ VNRTYQMHLG GKQREYEFFG LTPDTEFLGT IMICVPTWAK ESAPYMCRVK

image



Background

Alternative Names	Interleukin-22 receptor subunit alpha-1 (IL-22 receptor subunit alpha-1) (IL-22R-alpha-1) (IL-22RA1) (Cytokine receptor class-II member 9) (Cytokine receptor family 2 member 9) (CRF2-9) (ZcytoR11)
--------------------------	--

Note

Product Name: Recombinant Human IL22RA1 Protein
Catalog #: PEH5511



For Research Use Only , Not for Diagnostic Use