


Product Name: Recombinant Human PLXNA2 Protein
Catalog #: PEH5514



Summary

Name	PLXA2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	O75051
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	TEILTVSGPP EGGTRVTIHG VNLGLDFSEI AHHVQVAGVP CTPLPGEYII AEQIVCEMGH ALVGTTS GPV RLCIGECKPE FMTKSHQQYT FVNPSVLSLN PIRGPESGGT MVTITGHYLG AGSSVAVYLG NQTCEFYGRS MSEIVCVSPP SSNGLGPVPV SVSVDRAHVD SNLQFEYIDD PRVQRIEPEW SIASGHTPLT ITGFNLDVIQ EPRIRVKFNG KESVNVCKVV NTTTTCLAP SLTTDYRPGL DTVERPDEFG FVFNNVQSLL IYNDTKFIYY PNPTFELLSP TGVLDQKPGS PIILKGKNLC PPASGGAKLN YTVLIGETPC AVTVSETQLL CEPPNLTGQH KVMVHVGGMV

image



Background

Alternative Names	Plexin-A2 (Semaphorin receptor OCT)
--------------------------	-------------------------------------

Product Name: Recombinant Human PLXNA2 Protein
Catalog #: PEH5514



Note

For Research Use Only , Not for Diagnostic Use