

Product Name: Recombinant Human MME Protein
Catalog #: PEH5518



Summary

Name	NEP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P08473
Host	E.coli
Species	Human
Molecular Mass	45 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EKALAIKERI GYPDDIVSND NKLNNEYLEL NYKEDEYFEN IIQNLKFSQS KQLKKLREKV DKDEWISGAA VVNAFYSSGR NQIVFPAGIL QPPFFSAQQS NSLNYGGIGM VIGHEITHGF DDNGRNFNKG GDLVDWWTQQ SASNFKEQSQ CMVYQYGNFS WDLAGGQHLN GINTLGENIA DNGGLGQAYR AYQNYIKKNG EEKLLPGLDL NHKQLFFLNF AQVWCGTYRP EYAVNSIKTD VHSPGNFRRI GTLQNSAEFS EAFHCRKNSY MNPEKKCRVW

image



Background

Alternative Names	Neprilysin (EC 3.4.24.11) (Atriopeptidase) (Common acute lymphocytic leukemia antigen) (CALLA) (Enkephalinase) (Neutral endopeptidase 24.11) (NEP) (Neutral endopeptidase) (Skin fibroblast elastase) (SFE) (CD antigen)
--------------------------	--

Product Name: Recombinant Human MME Protein
Catalog #: PEH5518



CD10)

Note

For Research Use Only , Not for Diagnostic Use