



Product Name: Recombinant Human PLAGL1 Protein
Catalog #: PEH5519



Summary

Name	PLAL1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UM63
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	PYKCVQPD CG KAFVSR YKLM RHMATHSPQK SHQCAHCEKT FNRKDHLKNH LQTHDPNKMA FGCEECGKKY

image



Background

Alternative Names	Zinc finger protein PLAGL1 (Lost on transformation 1) (LOT-1) (Pleiomorphic adenoma-like protein 1) (Tumor suppressor ZAC)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use