


Product Name: Recombinant Human PRSS27 Protein
Catalog #: PEH5520



Summary

Name	PRS27_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9BQR3
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	MVGGQDTQEG EWPWQVSIQR NGSHFCGGS LIAEQWVLTA HCFRNTSETS LYQVLLGARQ LVQPGPHAMY ARVRQVESNP LYQGTASSAD VALVELEAPV PFTNYILPVC LPDPSVIFET GMNCWVTGWG SPSEEDLLPE PRILQKLAVP IIDTPKCNLL YSKDTEFGYQ PKTIKNDMLC AGFEEGKKDA CKGDSGGPLV CLVGQSWLQA GVISWGEGCA RQNRPGVYIR VTAHHNWIHR IIPKLQFQPA RLGGQK

image



Background

Alternative Names	Serine protease 27 (EC 3.4.21.-) (Marapsin) (Pancreasin)
--------------------------	--

Product Name: Recombinant Human PRSS27 Protein
Catalog #: PEH5520



Note

For Research Use Only , Not for Diagnostic Use