


Product Name: Recombinant Human STT3A Protein
Catalog #: PEH5526



Summary

Name	STT3A_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P46977
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LVDNNTWNNT HISRVGQAMA STEEKAYEIM RELDVSIVLV IFGGLTGYSS DDINKFLWMV RIGGSTDTGK HIKENDYYTP TGEFRVDREG SPVLLNCLMY KMCYYRFGQV YTEAKRPPGF DRVRNAEIGN KDFELDVLEE AYTTEHWLVR IYKVKDLDNR

image



Background

Alternative Names	Dolichyl-diphosphooligosaccharide--protein glycosyltransferase subunit STT3A (Oligosaccharyl transferase subunit STT3A) (STT3-A) (EC 2.4.99.18) (B5) (Integral membrane protein 1) (Transmembrane protein TMC)
--------------------------	--

Note

Product Name: Recombinant Human STT3A Protein
Catalog #: PEH5526



For Research Use Only , Not for Diagnostic Use