


Product Name: Recombinant Human POLR2J2 Protein
Catalog #: PEH5540



Summary

Name	RPB1B_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9GZM3
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKDTKVPKAC LFTINKEDHT LGNIIKSQLL KDPQVLFAGY KVPHPLEHKI IIRVQTTPDY SPQEAFTNAI TDLISELSLL

image



Background

Alternative Names	DNA-directed RNA polymerase II subunit RPB11-b1 (RNA polymerase II subunit B11-b1) (RPB11b1) (DNA-directed RNA polymerase II subunit J2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use