


Product Name: Recombinant Human POLR2J3 Protein
Catalog #: PEH5541



Summary

Name	RPB1C_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H1A7
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	NKDTKVPNAC LFTMNKEDHT LGNIIKSQLL KDPQVLFAGY KVPHPLEHKI IIRVQTTPDY SPQEAFTNAI TDLISELSLL

image



Background

Alternative Names	DNA-directed RNA polymerase II subunit RPB11-b2 (RNA polymerase II subunit B11-b2) (RPB11b2) (DNA-directed RNA polymerase II subunit J3)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use