


Product Name: Recombinant Human SLC38A3 Protein
Catalog #: PEH5555



Summary

Name	S38A3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q99624
Host	E.coli
Species	Human
Molecular Mass	9 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ELVPNGKHSE GLLPVITPMA GNQRVEDPAR SCMEGKSFLQ KSPSKEPHFT DFEGKTSFGM

image



Background

Alternative Names	Sodium-coupled neutral amino acid transporter 3 (N-system amino acid transporter 1) (Na ⁺)-coupled neutral amino acid transporter 3) (Solute carrier family 38 member 3) (System N amino acid transporter 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use