


Product Name: Recombinant Human SLURP2 Protein
Catalog #: PEH5559



Summary

Name	SLUR2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P0DP57
Host	E.coli
Species	Human
Molecular Mass	20 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IWCHQCTGFG GCSHGSRLR DSTHCVTTAT RVLSTEDLP LVTKMCHIGC PDIPSLGLGP YVSIACCQTS LCNHD

image



Background

Alternative Names	Secreted Ly-6/uPAR domain-containing protein 2 (Secreted LY6/PLAUR domain-containing protein 2) (Secreted Ly-6/uPAR-related protein 2) (SLURP-2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use