


Product Name: Recombinant Human TMCO4 Protein
Catalog #: PEH5561



Summary

Name	TMCO4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q5TGY1
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GKHLAHILLS RQGRRPVTI IGFSLGARVI YFCLQEMAE KDCQGIIEDV ILLGAPVEGE AKHWEPFRKV VSGRIINGYC RGDWLLSFVY RTSSVQLRVA GLQPVLLQDR RVENVDLTSV VSGHLDYAKQ MDAILKAVGI RTKPGWDEKG LLLAPGCLPS EEPRQAAAAA SSGETPHQVG QTQGPISGDT SKLAMSTDPS QAQVPVGLDQ SEGASLPAAA SPERPPICSH GMDPNPLGCP DCACKTQGPS

image



Background

Alternative Names	Transmembrane and coiled-coil domain-containing protein 4
--------------------------	---

Note

Product Name: Recombinant Human TCMO4 Protein
Catalog #: PEH5561



For Research Use Only , Not for Diagnostic Use