


Product Name: Recombinant Human TMED6 Protein
Catalog #: PEH5564



Summary

Name	TMED6_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8WW62
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GDQPLFRGAD RYDFAIMIPP GGTECFWQFA HQTGYFYFSY EVQRTVGMSH DRHVAATAHN PQGFLIDTSQ GVRGQINFST QETGFYQLCL SNQHNHFGSV QVYLNFGVFY EGPETDHKQK ERKQLNDTLD AIEDGTQKVQ

image



Background

Alternative Names	Transmembrane emp24 domain-containing protein 6 (p24 family protein gamma-5) (p24gamma5)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use