


Product Name: Recombinant Human TMPPE Protein
Catalog #: PEH5566



Summary

Name	TMPPE_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6ZT21
Host	E.coli
Species	Human
Molecular Mass	30 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KTVEVPIHQL PASMNNLKIV LLSDIHLGPT VGRTKMEMFV RMVNVLEPDI TVIVGDLSDS EASVLRATA PLGQLHSHLG AYFVTGNHEY YTSDVSNWFA LLESLHVQPL HNENVKISAT RAQRGGGGSG SGSEDEDWIC LAGVDDIEAD ILHYSGHGMD LDKALEGCSP DHTIILLAHQ PLAAKRALQA RPDINLILSG HTHAGQIFPL NVAAYLLNPF FAGLYQVAQA TFVYVSPGTA YYGIPMRLGS RAEITELILQ

image



Background

Alternative Names	Transmembrane protein with metallophosphoesterase domain (EC 3.1.-.-)
--------------------------	---

Product Name: Recombinant Human TMPPE Protein
Catalog #: PEH5566



Note

For Research Use Only , Not for Diagnostic Use