


Product Name: Recombinant Human TMPRSS2 Protein
Catalog #: PEH5567



Summary

Name	TMPS2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O15393
Host	E.coli
Species	Human
Molecular Mass	28 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKCSNSGIEC DSSGTCINPS NWCDGVSHCP GGEDENRCVR LYGPNFILQV YSSQRKSWHP VCQDDWNENY GRAACRDMGY KNNFYSSQGI VDDSGSTFSM KLNTSAGNVD IYKKLYHSDA CSSKAVVSLR

image



Background

Alternative Names	Transmembrane protease serine 2 (EC 3.4.21.-) (Serine protease 10) [Cleaved into: Transmembrane protease serine 2 non-catalytic chain; Transmembrane protease serine 2 catalytic chain]
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TMPRSS2 Protein
Catalog #: PEH5567

