


Product Name: Recombinant Human TNFSF8 Protein
Catalog #: PEH5587



Summary

Name	TNFL8_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P32971
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SPDNVPLKGG NCSDLLCIL KRAPFKKSWA YLQVAKHLNK TKLSWNKDCI LHGVRYQDGN LVIQFPGLYF IICQLQFLVQ CPNNSVDLKL ELLINKHIKK QALVTVCESG MQTKHVYQNL SQFLLDYLQV NTTISVNVDT FQYIDTSTFP LENVLSIFLY

image



Background

Alternative Names Tumor necrosis factor ligand superfamily member 8 (CD30 ligand) (CD30-L) (CD antigen CD153)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human TNFSF8 Protein
Catalog #: PEH5587

