

Product Name: Recombinant Human TM9SF4 Protein
Catalog #: PEH5588



Summary

Name	TM9S4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q92544
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EIKAVKLTSS RTQLPYEYYS LPFCQPSKIT YKAENLGEVL RGDRIVNTPF QVLMNSEKKC EVLCSQSNKP VTLTVEQSRL VAERITEDYY VHLIADNLPV ATRLELYSNR DSDDKKKEKD VQFEHGYRLG FTDVNKIYLH NHLSFILYYH REDMEDQEH TYRVVRFEEVI PQSIRLEDLK ADEKSSCTLP EGTNSSPQEI

image

SDS

Background

Alternative Names Transmembrane 9 superfamily member 4 (Tumor cannibalism associated protein 1)

Note

For Research Use Only , Not for Diagnostic Use