


Product Name: Recombinant Human SCGB2A1 Protein
Catalog #: PEH5593



Summary

Name	SG2A1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	O75556
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DSGCKLLEDM VEKTNSDIS IPEYKELLQE FIDSDAAAEA MGKFKQCFLN QSHRTLKNFG LMMHTVYDSI WCNMKS

image



Background

Alternative Names	Mammaglobin-B (Lacryglobin) (Lipophilin-C) (Mammaglobin-2) (Secretoglobulin family 2A member 1)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use