


Product Name: Recombinant Human SPOPL Protein
Catalog #: PEH5613



Summary

Name	SPOPL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6IQ16
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPECRLAEDL GNLWENTRFT DCSFFVRGQE FKAHKSVLAA RSPVFNAMFE HEMEESKKNR VEINDLDPEV FKEMMRFIYT GRAPNLDKMA DNLLAAADKY ALERLKVMCE EALCSNLSVE NVADTLVLAD LHSAEQLKAQ AIDFINRCSV LRQLGCKDGG NWNSNQATDI METSGWKSMI QSHPHLVAEA FRALASQCP QFGIPRKRLK

image



Background

Alternative Names Speckle-type POZ protein-like (HIB homolog 2) (Roadkill homolog 2)

Note

Product Name: Recombinant Human SPOPL Protein
Catalog #: PEH5613



For Research Use Only , Not for Diagnostic Use