


Product Name: Recombinant Human PRSS3 Protein
Catalog #: PEH5617



Summary

Name	TRY3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P35030
Host	E.coli
Species	Human
Molecular Mass	25 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IVGGYTCEEN SLPYQVSLNS GSHFCGGLI SEQWVVSAAH CYKTRIQVRL GEHNIKVLEG NEQFINAAKI IRHPKYNRDT LDNDIMLIK L SSPAVINARV STISLPTTPP AAGTECLISG WGNTLSFGAD YPDELKCLDA PVLQTAECKA SYPGKITNSM FCVGFLEGGK DSCQRDSGGP VVCNGQLQGV VSWG HGCAWK NRPGVYTKVY NYVDWIKDTI AANS

image



Background

Alternative Names	Trypsin-3 (EC 3.4.21.4) (Brain trypsinogen) (Mesotrypsin) (Mesotrypsinogen) (Serine protease 3) (Serine protease 4) (Trypsin III) (Trypsin IV)
--------------------------	--

Note

Product Name: Recombinant Human PRSS3 Protein
Catalog #: PEH5617



For Research Use Only , Not for Diagnostic Use