


Product Name: Recombinant Human LRRC4 Protein
Catalog #: PEH5627



Summary

Name	LRRC4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q9HBW1
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	DLNISEGRMA ELKCRTPPMS SVKWLLPNGT VLSHASRHRP ISVLNDGTLN FSHVLLSDTG VYTCMVTNVA GNSNASAYLN VSTAEINTSN YSFFTTVTVE TTEISPEDTT RKYKVPPTTS TGYQPAYTTS TTVLIQTTRV PKQVAVPATD TTDKMQTSLD

image



Background

Alternative Names	Leucine-rich repeat-containing protein 4 (Brain tumor-associated protein BAG) (Nasopharyngeal carcinoma-associated gene 14 protein) (Netrin-G2 ligand) (NGL-2)
--------------------------	--

Note

Product Name: Recombinant Human LRRC4 Protein
Catalog #: PEH5627



For Research Use Only , Not for Diagnostic Use