


Product Name: Recombinant Human LSAMP Protein
Catalog #: PEH5631



Summary

Name	LSAMP_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q13449
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SVDFNRGTDN ITVRQGDTAI LRCVVEDKNS KVAWLNRSI IFAGHDKWSL DPRVELEKRH SLEYSRLRIQK VDVEDEGSYT CSVQTQHEPK TSQVYLIVQV PPKISNISSD VTVNEGSNVT LVCMANGRPE PVITWRHLTP TGREFEGEEE YLEILGITRE QSGKYECKAA NEVSSADVQK VKVTVNYPPT ITESKSNEAT TGRQASKCE ASAVPAPDFE WYRDDTRINS ANGLEIKSTE GQSSLTVTNV TEEHYGNITC VAANKLGVTN

image



Background

Alternative Names Limbic system-associated membrane protein (LSAMP) (IgLON family member 3)

Product Name: Recombinant Human LSAMP Protein
Catalog #: PEH5631



Note

For Research Use Only , Not for Diagnostic Use