


Product Name: Recombinant Human MARCHF11 Protein
Catalog #: PEH5636



Summary

Name	MARHB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	A6NNE9
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	EPPPQPPPPP PTPPPGEPV PVPAAPRYLP PLPASPETPE RAAGPSEPLG EVAPRCRGAD ELPPPPLPLQ PAGQEVAAAG DSGEGPRRLP EAAAAGGGPG ESEAGAGGER ERRGAGDQPE TRSVCSSRSS SSGGGDQRAG HHHHHHHQPIC KICFQGAEQG ELLNPCRCDG

image



Background

Alternative Names	E3 ubiquitin-protein ligase MARCHF11 (EC 2.3.2.27) (Membrane-associated RING finger protein 11) (Membrane-associated RING-CH protein XI) (MARCH-XI) (RING-type E3 ubiquitin transferase MARCHF11)
--------------------------	---

Note



Product Name: Recombinant Human MARCHF11 Protein
Catalog #: PEH5636

For Research Use Only , Not for Diagnostic Use