


Product Name: Recombinant Human CHCHD1 Protein
Catalog #: PEH5637



Summary

Name	CHCH1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96BP2
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ARFGNPRKPV LKPNKPLILA NRVGERRREK GEATCITEMS VMMACWKQNE FRDDACRKEI QGFLDCAARA

image



Background

Alternative Names	Coiled-coil-helix-coiled-coil-helix domain-containing protein 1 (28S ribosomal protein S37, mitochondrial) (MRP-S37) (Mitochondrial small ribosomal subunit protein mS37) (Nuclear protein C2360)
--------------------------	---

Note

For Research Use Only , Not for Diagnostic Use