


Product Name: Recombinant Human ANOS1 Protein
Catalog #: PEH5643



Summary

Name	KALM_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P23352
Host	E.coli
Species	Human
Molecular Mass	26 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QRRRPTRPLE VGAPFYQDGQ LQVKVYWKKT EDPTVNRVHV RWFPEACAHN RTTGSEASSG MTHENYIILQ DLSFCKYKV TVQPIRPKSH SKAEAVFFTT PPCSALKGKS

image



Background

Alternative Names Anosmin-1 (Adhesion molecule-like X-linked) (Kallmann syndrome protein)

Note

For Research Use Only , Not for Diagnostic Use