



Summary

Name	HPDL_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q96IR7
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IAFHVPAGQP LARNLQRLFQ FQPLASREVD GWRQLALRSG DAVFLVNEGA GSGEPLYGLD PRHAVPSATN LCFDVADAGA ATRELAALGC SVPVPPVRVR DAQGAATYAV VSSPAGILSL TLLERAGYRG PFLPGFRPVS SAPGPGWVSR VDHLTLACTP GSSPTLLRWF HDCLGFCHLP LSPGEDPELG LEMTAGFGLG GLRLTALQAA PGSIVPTLVL AESLPGATTR QDQVEQFLAR HKGPGLQHVQ LYTPNIVEAT EGVATAGGQF LAPPGAYYQQ PGKERQIRAA GHEPHLLARQ GILLDGDGKGL FLLQVFTKSL FTEDTFFLEL IQRQGATGFG QGNIRALWQS

image



Background

Alternative Names	4-hydroxyphenylpyruvate dioxygenase-like protein (HPD-like protein) (EC 1.13.-.-) (Glyoxalase domain-containing protein 1)
--------------------------	--

Product Name: Recombinant Human HPDL Protein
Catalog #: PEH5651



Note

For Research Use Only , Not for Diagnostic Use