


Product Name: Recombinant Human MMP9 Protein
Catalog #: PEH5652



Summary

Name	MMP9_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P14780
Host	E.coli
Species	Human
Molecular Mass	21 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	VPDLGRFQTF EGDLEKWHHHN ITYWIQNYSE DLPRVIDDA FARAFALWSA VTPLTFTRVY SRDADIVIQF GVAEHGDGYP FDGKDGLLAH AFPPGPGIQG DAHFDDDELW SLGKGVVPT RFGNADGAAC HFPFIFEGRS YSACTTDGRS DGLPWCSTTA NYDTDDRFGF

image



Background

Alternative Names	Matrix metalloproteinase-9 (MMP-9) (EC 3.4.24.35) (92 kDa gelatinase) (92 kDa type IV collagenase) (Gelatinase B) (GELB) [Cleaved into: 67 kDa matrix metalloproteinase-9; 82 kDa matrix metalloproteinase-9]
--------------------------	---

Note

Product Name: Recombinant Human MMP9 Protein
Catalog #: PEH5652



For Research Use Only , Not for Diagnostic Use