


Summary

Name	ABRX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q6UWZ7
Host	E.coli
Species	Human
Molecular Mass	32 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SGFVLGALAF QHLNTDSDTE GFLLGEVKGE AKNSITDSQM DDVEVVYTID IQKYIPCYQL FSFYNSSGEV NEQALKKILS NVKKNVVGWY KFRRHSDQIM TFRERLLHKN LQEHFSNQDL VFLLLTPSII TESCSTHRLE HSLYKPKQGL FHRVPLVVAN LGMSEQLGYK TVSGSCMSTG FSRAVQTHSS KFFEEDGSLK EVHKINEMYA SLQEELKSIC KKVEDSEQAV DKLVKDVNRL KREIEKRRGA QIQAAREKNI

image



Background

Alternative Names	BRCA1-A complex subunit Abraxas 1 (Coiled-coil domain-containing protein 98) (Protein FAM175A)
--------------------------	--

Product Name: Recombinant Human ABRAXAS1 Protein
Catalog #: PEH5655



Note

For Research Use Only , Not for Diagnostic Use