


Product Name: Recombinant Human IL20RB Protein
Catalog #: PEH5662



Summary

Name	I20RB_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	Q6UXL0
Host	E.coli
Species	Human
Molecular Mass	23 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	HLVIELEDLG PQFEFLVAYW RREPGAEHV KMVRS GGIPV HLETMEPGAA YCVKAQTFVK AIGRYSAFSQ TECVEVQGEA

image



Background

Alternative Names Interleukin-20 receptor subunit beta (IL-20 receptor subunit beta) (IL-20R-beta) (IL-20RB) (Fibronectin type III domain containing 6) (FNDC6) (IL-20R2)

Note

For Research Use Only , Not for Diagnostic Use