


Product Name: Recombinant Human ITGB5 Protein
Catalog #: PEH5663



Summary

Name	ITB5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P18084
Host	E.coli
Species	Human
Molecular Mass	40 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	YYLMDLSLSM KDDLNIRSL GTKLAEEMRK LTSNFR LGFG SFVDKDISPF SYTAPRYQTN PCIGYKLFPN CVPSFGFRHL LPLTDRVDSF NEEVRKQRVS RNRDAPEGGF DAVLQAAVCK EKIGWRKDAL HLLVFTTDDV PHIALDGKLG GLVQPHDGQC HLNEANEYTA SNQMDYPSLA LLGEKLAENN INLIFAVTKN HYMLYKNFTA LIPGTTVEIL DGDSKNIIQL

image



Background

Alternative Names Integrin beta-5

Note

Product Name: Recombinant Human ITGB5 Protein
Catalog #: PEH5663



For Research Use Only , Not for Diagnostic Use