


Product Name: Recombinant Human LSM11 Protein
Catalog #: PEH5670



Summary

Name	LSM11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P83369
Host	E.coli
Species	Human
Molecular Mass	11 kDa
Buffer	Lyophilized from Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 12 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SKSAVEDSTL SRYSTSTWK LASVWGRADT GRGSHKRSRS VPSSLQASAR EESRSELSGR TTRTDGSSVG GTFSRATTLS RGQSRKKKRK

image



Background

Alternative Names U7 snRNA-associated Sm-like protein LSm11

Note

For Research Use Only , Not for Diagnostic Use