


Product Name: Recombinant Human DHX32 Protein
Catalog #: PEH5673



Summary

Name	DHX32_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q7L7V1
Host	E.coli
Species	Human
Molecular Mass	29 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	SSSEKRYFPE SLDSSDGDEE EVLACEDLEL NPFDGLPYSS RYKLLKERE DLPIWKEKYS FMENLLQNQI VIVSGDAKCG KSAQVPQWCA EYCLSIHYQH GGVICTQVHK QTVVQLALRV ADEMDVNIGH EVGYVIPFEN CCTNETILRY CTDDMLQREM MSNPFLGSYG VIILDDIHER SIATDVLLGL LKDVLLARPE LKLIINSSPH LISKLNSYYG NVPVIEVKNK

image



Background

Alternative Names	Putative pre-mRNA-splicing factor ATP-dependent RNA helicase DHX32 (EC 3.6.4.13) (DEAD/H box 32) (DEAD/H helicase-like protein 1) (DHLP1) (DEAH box protein 32) (HuDDX32)
--------------------------	---

Product Name: Recombinant Human DHX32 Protein
Catalog #: PEH5673



Note

For Research Use Only , Not for Diagnostic Use