

Product Name: Recombinant Human EPHB2 Protein
Catalog #: PEH5677



Summary

Name	EPHB2_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P29323
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	IMHQVSRTVD SITLSWSQPD QPNGVILDYE LQYYEKELSE YNATAIKSPT NTVTVQGLKA GAIYVFQVRA RTVAGYGRYS GKMYFQTMTE AEYQTSIQEK

image



Background

Alternative Names	Ephrin type-B receptor 2 (EC 2.7.10.1) (Developmentally-regulated Eph-related tyrosine kinase) (ELK-related tyrosine kinase) (EPH tyrosine kinase 3) (EPH-like kinase 5) (EK5) (hEK5) (Renal carcinoma antigen NY-REN-47) (Tyrosine-protein kinase TYRO5) (Tyrosine-protein kinase receptor EPH-3) [Cleaved into: EphB2/CTF1; EphB2/CTF2]
--------------------------	---

Note

Product Name: Recombinant Human EPHB2 Protein
Catalog #: PEH5677



For Research Use Only , Not for Diagnostic Use