


Product Name: Recombinant Human FMO4 Protein
Catalog #: PEH5686



Summary

Name	FMO4_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P31512
Host	E.coli
Species	Human
Molecular Mass	15 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	FERSDDIGGL WKFTESSKDG MTRVYKSLVT NVCKEMSCYS DFPFHEDYPN FMNHEKFDY LQEFAEHFDL LKYIQFKTTV CSITKRPDFS ETGQWDVVTE TEGKQNRAVF

image



Background

Alternative Names Dimethylaniline monooxygenase [N-oxide-forming] 4 (EC 1.14.13.8) (Dimethylaniline oxidase 4) (Hepatic flavin-containing monooxygenase 4) (FMO 4)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human FMO4 Protein
Catalog #: PEH5686

