


Summary

Name	LAMP5_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9UJQ1
Host	E.coli
Species	Human
Molecular Mass	24 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	QEVENLSGLS TNPEKDIFVV RENGTTCLMA EFAAKFIVPY DVWASNIVDL ITEQADIALT RGAEVKGRCG HSQSELQVFW VDRAVALKML FVKESHNMSK GPEATWRLSK VQFVYDSSEK THFKDAVSAG KHTANSHHLS ALVTPAGKSY ECQAQQTISL ASSDPQKTVT MILSAVHIQP FDIISDFVFS EEHKCPVDER

image



Background

Alternative Names	Lysosome-associated membrane glycoprotein 5 (Brain and dendritic cell-associated LAMP) (Brain-associated LAMP-like protein) (BAD-LAMP) (Lysosome-associated membrane protein 5) (LAMP-5)
--------------------------	--

Note

Product Name: Recombinant Human LAMP5 Protein
Catalog #: PEH5689



For Research Use Only , Not for Diagnostic Use