


**Product Name: Recombinant Human KDM5C Protein**  
**Catalog #: PEH5690**



## Summary

<b>Name</b>	KDM5C_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	P41229
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	14 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	ELMMEGDLLE VTLDENHSIW QLLQAGQPPD LERIRTLEL EKAERHGSRA RGRALERRRR RKVDRGGEGD DPAREELEPK RVRSSGPEAE EVQEEEELEE

## image



## Background

**Alternative Names** Lysine-specific demethylase 5C (EC 1.14.11.67) (Histone demethylase JARID1C) (Jumonji/ARID domain-containing protein 1C) (Protein SmcX) (Protein Xe169) ([histone H3]-trimethyl-L-lysine(4) demethylase 5C)

## Note

For Research Use Only , Not for Diagnostic Use