



Product Name: Recombinant Human LNX1 Protein
Catalog #: PEH5691



Summary

Name	LNX1_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q8TBB1
Host	E.coli
Species	Human
Molecular Mass	18 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LLVTCPFREH CTQVLQRCDL EHHFQTCKG ASHYGLTKDR KRRSQDGCPD GCASLTATAP SPEVSAATI SLMTDEPLD NPAYVSSAED GQPAISPVDS GRSNRTRARP FERSTIRSRS FKKINRALSV LRRTKSGSAV ANHADQGREN

image



Background

Alternative Names E3 ubiquitin-protein ligase LNX (EC 2.3.2.27) (Ligand of Numb-protein X 1) (Numb-binding protein 1) (PDZ domain-containing RING finger protein 2) (RING-type E3 ubiquitin transferase LNX)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human LNX1 Protein
Catalog #: PEH5691

