


**Product Name: Recombinant Human MMP17 Protein**  
**Catalog #: PEH5701**



## Summary

<b>Name</b>	MMP17_HUMAN
<b>Purity</b>	Greater than 85% as determined by reducing SDS-PAGE
<b>Endotoxin level</b>	Please contact with the lab for this information
<b>Tag</b>	C-His
<b>Accession #</b>	Q9ULZ9
<b>Host</b>	E.coli
<b>Species</b>	Human
<b>Molecular Mass</b>	21 kDa
<b>Buffer</b>	Phosphate buffered saline
<b>Shipping</b>	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
<b>Stability&amp;Storage</b>	Store at $\leq -70^{\circ}\text{C}$ , stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$ , stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
<b>Sequence</b>	RTSDHKIVFF KGDYRWVFKD NNVEEGYPRP VSDFSLPPGG IDAAFSWAHN DRTYFFKDQL YWRYDDHTRH MDPGYPAQSP LWRGVPSTLD DAMRWSDGAS YFFRGQEYWK VLDGELEVAP GYPQSTARDW LVCGDSQADG SVAAGVDAAE GPRAPPGQHD QSRSEGGYEV

## image



## Background

<b>Alternative Names</b>	Matrix metalloproteinase-17 (MMP-17) (EC 3.4.24.-) (Membrane-type matrix metalloproteinase 4) (MT-MMP 4) (MTMMP4) (Membrane-type-4 matrix metalloproteinase) (MT4-MMP) (MT4MMP)
--------------------------	---

## Note

**Product Name: Recombinant Human MMP17 Protein**  
**Catalog #: PEH5701**



---

For Research Use Only , Not for Diagnostic Use