

Product Name: Recombinant Human DDX50 Protein
Catalog #: PEH5703



Summary

Name	DDX50_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9BQ39
Host	E.coli
Species	Human
Molecular Mass	41 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	ESQKKERQKS DRRKSRHHYD SDEKSETREN GVTDDLDA PK AKKSKMKEKL NGDTEEGFNR LSDEFKSHK SRRKDLPNGD IDEYEKSKR VSSLDSTHK SSDNKLEETL TREQKEGAFS NFPISEETIK LLKGRGVTYL FPIQVKTFGP VYEGKDLIAQ ARTGTGKTFS FAIPLIERLQ RNQETIKSR SPKVLVLAPT RELANQVAKD FKDITRKLSV ACFYGGTSYQ SQINHIRNGI DILVGTPGRI KDHLLQSGRLD LSKLRHVVDL EVDQMLDLGF AEQVEDIHE SYKTDSEDP QTLLFSATCP QWVYKVAKKY MKSRYEQVDL VGKMTQKAAT

image

SDS

Background

Alternative Names	ATP-dependent RNA helicase DDX50 (EC 3.6.4.13) (DEAD box protein 50) (Gu-beta) (Nucleolar protein Gu2)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use