


Product Name: Recombinant Human EPHB3 Protein
Catalog #: PEH5705



Summary

Name	EPHB3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	P54753
Host	E.coli
Species	Human
Molecular Mass	13 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LHSSGSSLT LSWAPPERPN GVILDYEMKY FEKSEGIAS T VTSQMNSVQL DGLRPDARYV VQVRARTVAG YGQYSRPAEF ETTSERGSGA

image



Background

Alternative Names Ephrin type-B receptor 3 (EC 2.7.10.1) (EPH-like tyrosine kinase 2) (EPH-like kinase 2) (Embryonic kinase 2) (EK2) (hEK2) (Tyrosine-protein kinase TYRO6)

Note

For Research Use Only , Not for Diagnostic Use