


Product Name: Recombinant Human FUT11 Protein
Catalog #: PEH5706



Summary

Name	FUT11_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q495W5
Host	E.coli
Species	Human
Molecular Mass	19 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	LAYKQPGGIT NQFLDSLKH REWGVNDPLL PNYLNGFECF VCDYELARLD AEKAHAASPG DSPVFEPHIA QPSHMDCPVP TPGFGNVEEI PENDSWKEMW LQDYWQGLDQ GEALTAMIHN NETEQTKFWD YLHEIFMKRQ

image



Background

Alternative Names	Alpha-(1,3)-fucosyltransferase 11 (EC 2.4.1.-) (Fucosyltransferase XI) (Fuc-TXI) (FucT-XI) (Galactoside 3-L-fucosyltransferase 11) (Fucosyltransferase 11)
--------------------------	--

Note

For Research Use Only , Not for Diagnostic Use