


Product Name: Recombinant Human HTRA3 Protein
Catalog #: PEH5707



Summary

Name	HTRA3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	N-SUMO & C-His
Accession #	P83110
Host	E.coli
Species	Human
Molecular Mass	31 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs.Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt.Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	GRELGLRDS D MDYIQTDAII NYGNSGGPLV NLDGEVIGIN TLKVTAGISF AIPSDRITRF LTEFQDKQIK DWKKRFIGIR MRTITPSLVD ELKASNPDFP EVSSGIYVQE VAPNSPSQRG GIQDGDIIVK VNGRPLVDSS ELQEAVLTES PLLLEVRRGN

image



Background

Alternative Names Serine protease HTRA3 (EC 3.4.21.-) (High-temperature requirement factor A3) (Pregnancy-related serine protease)

Note

For Research Use Only , Not for Diagnostic Use

Product Name: Recombinant Human HTRA3 Protein
Catalog #: PEH5707

