


Product Name: Recombinant Human MINDY3 Protein
Catalog #: PEH5715



Summary

Name	MINY3_HUMAN
Purity	Greater than 85% as determined by reducing SDS-PAGE
Endotoxin level	Please contact with the lab for this information
Tag	C-His
Accession #	Q9H8M7
Host	E.coli
Species	Human
Molecular Mass	10 kDa
Buffer	Phosphate buffered saline
Shipping	The product is shipped on dry ice/polar packs. Upon receipt, store it immediately at the temperature listed below.
Stability&Storage	Store at $\leq -70^{\circ}\text{C}$, stable for 6 months after receipt. Store at $\leq -70^{\circ}\text{C}$, stable for 3 months under sterile conditions after opening. Please minimize freeze-thaw cycles.
Sequence	KSSWRDCSEE EQKELLCHTL CDILESACCD HSGSYCLVSW LRGKTTEETA SISGSPAESS CQVEHSSALA

image



Background

Alternative Names Ubiquitin carboxyl-terminal hydrolase MINDY-3 (EC 3.4.19.12) (Dermal papilla-derived protein 5) (Deubiquitinating enzyme MINDY-3) (Protein CARP)

Note

For Research Use Only , Not for Diagnostic Use